
## MR Safe Eye Tracker





## LiveTrack AV for fMRI

#### Affordable 60Hz video eye tracking.

LiveTrack delivers real time estimates of eye rotation, direction of gaze coordinates and pupil size. Our close to the eye imaging arrangement is compact, easy to set up and delivers more robust tracking than a remote system.

Monocular and binocular configurations are available for a range of different head coils.

LiveTrack integrates with our BOLDscreen MR Safe displays, audio system, and the other *Made for fMRI* devices from our range.

Designed for fMRI by CRS

www.crsltd.com/livetrack-fmri

# LiveTrack AV for fMRI

#### Dedicated eye tracking hardware

Bespoke embedded hardware identifies the pupil and corneal reflection, then calculates eye rotation. Video images, tracked eye position and pupil size data relayed to host computer via USB interface.

#### Monitor participant's behaviour via live video

Live video stream from each eye tracking camera is produced concurrently with the data. Get instant feedback to know where your participant is looking.

#### Very easy integration

Works with Presentation, ePrime, MATLAB, Python and your own custom software. Compatible with Windows, Mac OS X and GNU/Linux.

### Cold mirror optical system

LiveTrack

MR safe cameras, with integrated infrared illumination sources, are positioned behind a large MR safe cold mirror. Cold mirrors are special optical components that are coated with a dichroic filter that reflects visible light while transmitting infrared wavelengths.

This allows the cameras to be optimally aligned in direct sight of the participant's eyes, but at the same time hidden behind the mirror outside the participant's view. The cameras can unobtrusively image each eye using infrared light through the mirror.

The participant sees the visible light image reflected from the cold mirror, and can therefore view a display placed outside the magnet. Our BOLDscreen MR Safe LCD displays are ideal for this purpose: they can be placed right at the opening of the bore, to give maximum field of view.

#### Fast setup, monocular and binocular tracking

Mounting the cameras and illumination directly on the head coil avoids all the difficulties of remote camera systems.

A range of opto-mechanical mounts are available to support different head coils; monocular and binocular configurations for tracking with one or two cameras.

Cambridge Research Systems Tel: +44 1634 720707 USA/Canada Toll Free: 1 866 846 2929 Email: enquiries@crsltd.com

www.crsltd.com

### For more details: www.crsltd.com/livetrack-fmri





M0203

