



## EcoFlexMEA36

Temperature compatibility Dimensions (W x D x H)

Thickness of the electrode field

Track material and contact pads

Flexible Microelectrode Array with 36 Electrodes for Use with ME2100-HS32 Headstage or 32-Channel Miniature Preamplifier MPA32I for in vivo or in vitro applications.

Layout

Base material

Electrode diameter

(center to center)

Electrode height

Electrode material

Isolation material

Electrode impedance

Electrode layout grid

Ground electrodes

Software

MC\_Rack

Channel map

Number of recording electrodes

Number of reference electrodes

Multi Channel Experimenter

Interelectrode distance

Diameters of the holes

Weight

#### **Technical Specifications**

0 - 125 °C 37 mm x 30 mm 50 µm Polyimide (Kapton) < 10 g Gold (Au) 50 µm 300 µm 100 µm Planar Gold (Au) Polyimide (Kapton)  $< 150 \text{ k}\Omega$ 6 x 6 32 2 internal reference electrode (iR) 2 ground electrodes

**MEA** Configuration 1 dimensional or Configuration EcoFlexMEA36.cmp

## Cleaning

Rinse with distilled water, optional with ethanol 70 %. Do not autoclave EcoFlexMEAs made from Polyimide (Kapton).

Multi Channel Systems MCS GmbH Aspenhaustrasse 21 72770 Reutlingen Germany

Phone +49-7121-909 25- 0 +49-7121-909 25-11

sales@multichannelsystems.com www.multichannelsystems.com

© 2019 Multi Channel Systems MCS GmbH a division of Harvard Bioscience, Inc.

Product information is subject to change without notice

for grounding the bath.

mm

Connect the EcoFlexMEA36 directly to a ME2100-HS32 or 32-Channel Miniature

the MPA32I with the electrode field up.

The additional connector can be used for connecting a silver pellet or a silver wire

Preamplifier. Insert the EcoFlexMEA36 into

30 mm

6 mm

28 mm

mm

14 mm

37 mm

#### **Advantages**

- EcoFlexMEAs made of flexible polyimide (Kapton) are perfect for in vivo and specific in vitro applications.
- EcoFlexMEAs are very cost efficient and more robust than FlexMEAs from polyimide foil.
- The electrodes and tracks and contact pads are made of pure gold.

Fax





# **EcoFlexMEA36**

Flexible Microelectrode Array with 36 Electrodes for Use with ME2100-HS32 Headstage or 32-Channel Miniature Preamplifier MPA32I for *in vivo* or *in vitro* applications.

Layout



**Electrode field:** 1800 x 1800 μm,

32 recording electrodes 2 reference electrodes, 2 ground electrodes.

GND 1 is a large ground electrode connected to pin 1 of the MPA32I input connector. GND 2 is a second ground electrode connected to pin 36. The REF electrodes are reference electrodes connected to pin 2 and 35, respectively. Both ground inputs and both reference electrode inputs are equal, that is, they are connected to each other inside the standard MPA32I. Please see the MPA32I manual for details.



The numbers in the electrodes are the recording channel numbers that refer to the channel numbers of the data acquisition program. For MC\_Rack, please make sure that you have selected "Configuration" in the "Channel Layout" under "Data Source Setup" with a total number of 64 channels. In "Amplifier", please choose FA32I/S or FA64I/S and in "MEA" EcoFlex-MEA36. In Layout tab of the display, please click "Default Map".

### Pin Layout

EcoFlexMEA36 Input Pins

**Important:** The data acquisition channel map is constructed by looking on the back side of the electrodes, because the FlexMEA

electrodes are placed on the pre-

22 D5 20 E5

300 µm

GN

27 B6 25 C6 27 D6

300 µm

paration upside down!

5 B1 7 C1 9 D1 1



EcoFlexMEA36 Input Pins

Pin 1GND (Ground)Pin 2Reference inputPin 3 to 34Recording channels 1 to 32Pin 35ReferencePin 36GND (Ground)

Warning: The device may only be used together with the ME2100-HS32 headstage or th MPA32I (-Flex) from Multi Channel Systems MCS GmbH, and only for the specified purpose. Damage of the device and even injuries can result from improper use.

Multi Channel Systems MCS GmbH Aspenhaustrasse 21 72770 Reutlingen Germany Phone +49-7121-909 25- 0 Fax +49-7121-909 25-11

sales@multichannelsystems.com www.multichannelsystems.com © 2019 Multi Channel Systems MCS GmbH a division of Harvard Bioscience, Inc.

Product information is subject to change without notice.