




# FlexMEA72

Flexible Microelectrode Array with 72 electrodes for use via ADPT-FM-72 adapter with two 32-Channel Miniature Preamplifier MPA32I for in vivo and in vitro applications.

Layout



# Temperature compatibility Dimensions (W x D x H) Thickness of the electrode field Base material Weight Track material and contact pads Electrode diameter Interelectrode distance (center to center) Diameters of the holes Electrode height Electrode material Isolation material Electrode impedance Electrode layout grid Number of recording electrodes Number of reference electrodes Ground electrodes Software Multi Channel Experimenter MC\_Rack Channel map

#### **Technical Specifications**

- 0 125 °C 42 mm x 12 mm 12 µm Polyimide 2611 foil < 1 g Gold (Au) 100 µm 625 to 750 µm
- 100 µm Planar TiN electrodes (Titanium nitride) Polyimide 2611 foil  $< 150 \text{ k}\Omega$ 8 x 9 64 4 internal reference electrode (iR) 4 ground electrodes

**MEA** Configuration 1 dimensional or Configuration FlexMEA72.cmp

## **Advantages**

- FlexMEAs made of flexible polyimide 2611 foil, perfect for in vivo and specific in vitro applications, for example, whole-heart preparations.
- The titanium nitride electrodes have a diameter of 100 µm, and the distance between the electrodes is 625 to 750 µm.
- The polyimide foil is perforated with holes of 100 µm diameter, ensuring optimal tissue contact.

## **Cleaning and Sterilization**

Rinse with distilled water, optional with ethanol 70 %. Please do not use an ultrasonic bath! When using an autoclave at 121 °C, please make sure that the FlexMEA will not be exposed to the moisture. The FlexMEA itself must be dry and additionally sealed in a sterile package.


Multi Channel Systems MCS GmbH Aspenhaustrasse 21 72770 Reutlingen Germany

Phone +49-7121-909 25-0 +49-7121-909 25-11 Fax

sales@multichannelsystems.com www.multichannelsystems.com

© 2019 Multi Channel Systems MCS GmbH a division of Harvard Bioscience, Inc.


Product information is subject to change without notice



The numbers in the electrodes are the recording channel numbers that refer to the channel numbers in the MC\_Rack program. Please make sure that you have selected the linear (1 dimensional) channel layout or Configuration under "Data Source Setup" with a total number of 64 channels. Deselect the check box "Digital input channel" if you do not need it, otherwise one recording channel is missing! See the MC\_Rack manual or help for details. The letter digit code below is the electrode identifier and refers to the position of the electrode in the grid.

**Important:** The MC\_Rack channel map is constructed by looking on the back side of the electrodes, because the FlexMEA72 electrodes are placed on the preparation upside down!

Please make sure to connect the left hand side of the Flex-MEA72 contact pads (electrode side on top) via ADPT-FM-72 adapter to the first MPA32I (channels 33 to 64), and the right hand side to the second MPA32I (channels 1 to 32). Read the SC2x32 data sheet when using a signal collector.



MPA32I Channel 33 - 64 MPA32I Channel 1 - 32



Warning: The device may only be used together with the MPA32I (-Flex) from Multi Channel Systems MCS GmbH, and only for the specified purpose. Damage of the device and even injuries can result from improper use.

Multi Channel Systems MCS GmbH Aspenhaustrasse 21 72770 Reutlingen Germany Phone +49-7121-909 25- 0 Fax +49-7121-909 25-11 © 2019 Multi Channel Systems MCS GmbH a division of Harvard Bioscience, Inc.

sales@multichannelsystems.com www.multichannelsystems.com Product information is subject to change without notice.