



a division of Harvard Bioscience, Inc.


2200

Dimensions in µm

FlexMEA36

Flexible Microelectrode Array with 36 Electrodes for Use with ME2100-HS32 Headstage or 32-Channel Miniature Preamplifier via Adapter ADPT-FM-32.

Layout



Important: The channel map is constructed by looking on the back side of the electrodes, because the FlexMEA electrodes are placed on the preparation upside down!

Technical Specifications

Temperature compatibility Dimensions (W x D x H)

Thickness of the electrode field

Base material

Weight

Track material and contact pads

Electrode diameter

Interelectrode distance

(center to center)

Diameters of the holes Electrode height

Electrode material

Isolation material

Electrode impedance

Electrode layout grid

Number of recording electrodes

Number of reference electrodes

Ground electrodes

Software

Multi Channel Experimenter

MC Rack Channel map 0 - 125 °C

41 mm x 11 mm

12 µm

Polyimide 2611 foil

< 1 q

Gold (Au)

30 µm

300 µm

50 µm

Planar

TiN electrodes (Titanium nitride)

Polyimide 2611 foil

< 150 $k\Omega$

6 x 6

2 internal reference electrode (iR)

2 ground electrodes

MEA Configuration

1 dimensional or Configuration

FlexMEA36.cmp

Advantages

The FlexMEA36 is made of flexible polyimide 2611 foil, perfect for in vivo and specific in vitro applications, for example, whole-heart preparations.

50d

- The titanium nitride electrodes have a diameter of 30 µm, and the distance between the electrodes is 300 µm.
- The polyimide foil is perforated with holes of 30 µm diameter, ensuring optimal tissue contact.

Cleaning and Sterilization

Rinse with distilled water, optional with ethanol 70 %. Please do not use an ultrasonic bath! When using an autoclave at 121 °C, please make sure that the FlexMEA36 will not be exposed to the moisture. The FlexMEA36 itself must be dry and additionally sealed in a sterile package.

Multi Channel Systems MCS GmbH Aspenhaustrasse 21 72770 Reutlingen Germany

Phone Fax

+49-7121-909 25- 0 +49-7121-909 25-11

sales@multichannelsystems.com www.multichannelsystems.com

© 2019 Multi Channel Systems MCS GmbH a division of Harvard Bioscience, Inc.

Product information is subject to change without notice






a division of Harvard Bioscience, Inc.




Flexible Microelectrode Array with 36 Electrodes for Use with ME2100-HS32 Headstage or 32-Channel Miniature Preamplifier via Adapter ADPT-FM-32.

Layout



Setup

The labelled side with NMI / A1 / A6 / F1 / F6 is the correct side with the contact pads and electrodes. It might be a bit confusing that the pads look stronger from the wrong side, but if you hold the FlexMEA36 into the light, you see that the pads have a 3-dimensional appearance only from the correct side. Please insert the FlexMEA upside down into the ADPT-FM-32 and connect the adapter to the MPA32I or to the ME2100-HS32



If you use more than one MPA32I and a MEA64-System, the signal collector SC2x32 leads the output channels of the second amplifier to channel number 33 - 64. Please see datasheet SC2x32 for details.

The numbers in the electrodes are the recording channel numbers that refer to the channel numbers in the Multi Channel Experimenter or MC_Rack program. The letter digit code below is the electrode identifier and refers to the position of the electrode in the grid. Please consider the orientation! For MC_Rack, please make sure that you have selected the linear (1 dimensional) channel layout under "Data Source Setup" with a total number of 32 channels. Deselect the check box "Digital input channel" if you do not need it, otherwise one recording channel is missing! See the MC_Rack or Multi Channel Experimenter manual or help for details.



Warning: The device may only be used together with the ME2100-HS32 or the MPA32I (-Flex) from Multi Channel Systems MCS GmbH, and only for the specified purpose. Damage of the device and even injuries can result from improper use.