



a division of Harvard Bioscience, Inc.

FlexMEA36-OM

Flexible Microelectrode Array with 36 Electrodes and Omnetics Connector for Use with the ME2100-HS32 Headstage or the 32-Channel Micro Preamplifier µPA32

Layout

Cleaning

Rinse with distilled water, optional with ethanol 70 %. Electrode field on this side

Sterilization: Do not autoclave or sterilize FlexMEAs by heat. These MEA types are not heat-stable and will be irreversibly damaged! Please do not use ultrasonic bath for cleaning.

Technical Specifications

10 - 50 °C Temperature compatibility

Dimensions (W x D x H) 45 mm x 13.8 mm x 1.8 mm

Thickness of the electrode field 12 µm

Base material Polyimide 2611 foil

Weight < 1 gTrack material and contact pads Gold (Au) Flectrode diameter 30 µm

Interelectrode distance (center to center)

Diameter of the holes 50 µm Electrode height Planar

Electrode material TiN electrodes (Titanium nitride)

300 µm

Isolation material Polyimide 2611 foil

Electrode impedance $< 150 \text{ k}\Omega$

6 x 6 Electrode layout grid Number of recording electrodes 32 Number of reference electrodes 2

Number of ground electrodes 2

Software

Multi Channel Experimenter linear

Advantages

- The FlexMEA36-OM is made of flexible polyimide 2611 foil, perfect for in vivo and specific in vitro applications, for example, whole-heart preparations.
- The titanium nitride electrodes have a diameter of 30 µm, and the distance between the electrodes is 300 µm.
- The polyimide foil is perforated with holes of 30 µm diameter, ensuring optimal tissue contact.



Warning: The FlexMEA36-OM may only be used together with the devices from Multi Channel Systems MCS GmbH, and only for the specified purpose. Damage of the device and even injuries can result from improper use.






a division of Harvard Bioscience, Inc.

FlexMEA36-OM

Flexible Microelectrode Array with 36 Electrodes and Omnetics Connector for Use with the ME2100-HS32 Headstage or the 32-Channel Micro Preamplifier µPA32


Layout

Important: The channel map is constructed by looking on the back side of the electrodes, because the FlexMEA36-OM electrodes are placed on the preparation upside down!



Numbering

Electrode Layout (electrodes facing down, seen from behind)



▼ Direction to the connector. The letter digit code below the electrode refers to the position of the electrode in the grid. The numbers inside are the linear channel numbers in Multi Channel Experimenter and MC_Rack.

Multiple Headstages

MC Rack


If you use more than one µPA32, the signal collector SC2x32 leads the output channels of the second amplifier to channel number 33 - 64. Please see datasheet SC2x32 for details.

Multi Channel Experimenter

If you use a ME2100-System with more than one headstage, the second/third/ fourth headstage will be displayed as channels B1 - B32 / C1 - C32 / D1 - D32.

Connector Layout

Labeled side of the connector.



Connector Specification

Manufacturer: Omnetics
Manufacturer Part No: A79022-001
Type: NDP-36-DD-GS

Connector with 36 pos, dual row, male, provided with 4 holes for guide posts.

 \bigotimes

Guide post

January 2020

Multi Channel Systems MCS GmbH Aspenhaustrasse 21 72770 Reutlingen Germany Phone +49-7121-909 25- 0 Fax +49-7121-909 25-11

sales@multichannelsystems.com www.multichannelsystems.com © 2020 Multi Channel Systems MCS GmbH a division of Harvard Bioscience, Inc.

Product information is subject to change without notice.