



# **CMOS-MEA Chip**

for Use with CMOS-MEA5000-System

**Technical Specifications** 



**Recording Sites** 

The CMOS-MEA chip is based on complementary metal oxide semiconductor (CMOS) technology

# 55 mm

55 mm

## Important!

CMOS-MEAs are not symmetrical! Please take care for the correct orientation of the chip. The round edge of the CMOS-MEA has to be in the front on the left side when looking directly to the open CMOS-MEA headstage.

#### **Sterilization**

Before using the CMOS-MEA chip, please sterilize the surface with UV radiation, for example, in a conventional flow box.

Warning: Autoclaving of CMOS-MEAs is possible, but not recommended as standard procedure.

Warning: It is not recommended to treat the CMOS chip in plasma cleaner! Extensive plasma treatment might

damage the CMOS-MEAs. Keep treatment short and low power: 60 s at 0.2 mbar and 5 to 10 W.

#### **Technical Specifications**

Temperature Dimensions (W x D x H)

Base material CMOS chip Track material Contact pads

Electrode diameter Interelectrode distance from center to center Active area

Flat area (around active area)

Recording electrodes

Stimulation electrodes

0 - 125 `°C 55 x 55 x 2 mm

PCB glass like surface Bonding wires and Au Au

8 µm

16 or 32 µm 1.04 mm x 1.04 mm or 2.08 mm x 2.08 mm rectangle: 8 mm x 8.5 mm or round: diameter 6 mm

4225 in 65 x 65 layout grid 1024 in 32 x 32 layout grid

Multi Channel Systems MCS GmbH Aspenhaustrasse 21 72770 Reutlingen Germany

Phone +49-7121-909 25-0 +49-7121-909 25-11 Fax

sales@multichannelsystems.com www.multichannelsystems.com

© 2021 Multi Channel Systems MCS GmbH a division of Harvard Bioscience, Inc.

May 2020

Product information is subject to change without notice.



on the CMOS-MEA chip. It has a flat rectangle area for placing the slice (8 mm x 8.5 mm) with influx and drain for intensive perfusion. Please use an external Ag/AgCl electrode to ground the bath.

Use the Culture Chamber CC with lid for cell cultures on the CMOS-MEA chip. It has a flat round area (diameter 6 mm) for the cell culture. It is possible to cultivated the cells on the chip in an incubator. A ground electrode is integrated in the CC, an external Ag/AgCl electrode is not necessary.

Please note: The lid has to be sealed with a semipermeable membran. We recommended ALA MEA-MEM-Sheets from ALA Scientific, Inc., which need to be ordered separately.



## Cleaning of the CMOS-MEA Chip

CMOS chips are gently cleaned with detergent Tickopur R36 (5%, Stamm/Berlin, 80 degC), and rinsed with ultrapure water (resistivity:18 MVcm). Fill hot Tickopur solution (80 °C) with a glass pipette into the CMOS MEA chamber onto the electrode field. Remove the Tickopur solution after 2 minutes by rinsing the chip with ultrapure water.

Use a cotton swab for carefully cleaning the surface, if necessary. Please do not damage the surface mechanically, otherwise the chip will be destroyed.

The following cleaning protocol was used in experiments with retina cells. Please read the paper which is online free available: Axonal Transmission in the Retina Introduces a Small Dispersion of Relative Timing in the Ganglion Cell Population Response from Guenther Zeck, Armin Lambacher, Peter Fromherz (2011).

Multi Channel Systems MCS GmbH Aspenhaustrasse 21 72770 Reutlingen Germany Phone+49-7121-909 25- 0Fax+49-7121-909 25-11

May 2020 © 2021 Multi Channel Systems MCS GmbH

sales@multichannelsystems.com www.multichannelsystems.com Product information is subject to change without notice.

a division of Harvard Bioscience, Inc.